



# Recombinant *Suillus bovinus* Tubulin beta-1 chain (TUBB1)

<b>Product Code</b>	CSB-MP728617SFAP
<b>Abbreviation</b>	TUBB1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6EVK8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Suillus bovinus</i> (Jersey cow bolete) ( <i>Boletus bovinus</i> )
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MREIVHIQTG QCGNQIGAKF WEVVSDEHGI ERDGLYKGTN DMQLERISVY YNEIGSNKYV PRAVLVDLEP GTMDSVRS GP LGGLFRPDNF IFGQSGAGNN WAKGHYTEGA ELVDSVLDVV RKEAEGTDCL QGFQITHSLG GGTGAGMGTL LMSKIREEYP DRMMCTYSVV PSPAVSDTVV EPYNATLSVH QLVENSDETFCIDNEALYDI CFRTLKLSTP TYGDLNHLVS FVMSGITTCL RFPGQLNSDL RKLAVNMVPF PRLHFFMTGF APLTARGSQQ YRAVTVPELT QQMFDKNNMM AASDPRHG RY LTVA AVFRGK VSMKEVEEQM QNVQNKNSAY FVEWIPNNVL SAQCDIPPRG VKMAVTFLGN STAIQELFKR VSDHFTAMFK RKAFLHWYTQ EGMDEMEFTE AESNMQDLIA EYQQYQDATV EEEAEYEEEEV PADEES
<b>Source</b>	Mammalian cell
<b>Target Names</b>	TUBB1
<b>Protein Names</b>	Recommended name: Tubulin beta-1 chain Alternative name(s): Beta-1-tubulin
<b>Expression Region</b>	1-446
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.