



Recombinant *Mesoplasma florum* Ribonuclease P protein component (rnpA)

Product Code	CSB-EP728634MEO
Abbreviation	rnpA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6F0D5
Product Type	Recombinant Protein
Immunogen Species	<i>Mesoplasma florum</i> (strain ATCC 33453 / NBRC 100688 / NCTC 11704 / L1) (<i>Acholeplasma florum</i>)
Purity	≥85% (SDS-PAGE)
Sequence	MKNKKIKN FEFQEISKQ EFHRNSAFVI YYSKNDKGYF RYGISVGKKL GNAVTRNKIK RQIRMMIQDQ IKILPEFSYD IVIARNRMM QNSFDQNKQE LNKLVVRFK
Source	<i>E.coli</i>
Target Names	rnpA
Protein Names	Recommended name: Ribonuclease P protein component Short name= RNase P protein Short name= RNaseP protein EC= 3.1.26.5 Alternative name(s): Protein C5
Expression Region	1-110
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.