



Recombinant Salmonella paratyphi A Allantoinase (allB)

Product Code	CSB-EP728664SBM-B
Abbreviation	allB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5PCG0
Product Type	Recombinant Protein
Immunogen Species	Salmonella paratyphi A (strain ATCC 9150 / SARB42)
Purity	≥85% (SDS-PAGE)
Sequence	MSFDLIKNG TVILENEARV IDIAVQGGKI AAIGENLGEA KNVLDTGLI VSPGMVDAHT HISEPGRTHW EGYETGTRAA AKGGITTMIE MPLNQLPATV DRETIELKFD AAKGKLTIDA AQLGGLVSYN LDRLHELDEV GVVGFKCFVA TCGDRGIDND FRDVNDWQFY KGAQKLGE MD QTVLVHCENA LICDELGEEA KREGRTAHD YVASRPVFTE VEAIRRVLYL AKAAGCRLHV CHISSPEGVE EVTRARQEGQ DVTCESEPHY FVLDTDQFEE IGTAKCSPP IRDQENQKGM WEKLFNGEID CLVSDHSPCP PEMKAGNIMQ AWGGIAGLQN CMDVMFDEAV QKRGMSLPMF GKLMATNAAD IFGLKHKGRI APGKDADLVF IQPDSSYVLK NEDLEYRHKV SPYVGR TIGA RITKTILRGD VIYDIEHGFP VPPKGQFILK HQQ
Source	E.coli
Target Names	allB
Protein Names	Recommended name: Allantoinase EC= 3.5.2.5 Alternative name(s): Allantoin-utilizing enzyme
Expression Region	1-453
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.