



Recombinant Salmonella paratyphi A Twin-arginine leader-binding protein DmsD (dmsD)

Product Code	CSB-MP728711SBM
Abbreviation	dmsD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5PHH8
Product Type	Recombinant Protein
Immunogen Species	Salmonella paratyphi A (strain ATCC 9150 / SARB42)
Purity	≥85% (SDS-PAGE)
Sequence	MTTFLQRDDF AVTARVLGAL FYSPESHET APLVQALLND DWQAQWPLDA EALAPVAVMF KTHSEESLPQ AWQRLFIGPY ALPSPWGSV WLDRESVLFG DSTLALRQWM RENGIQFEMQ QNEPEDHFGS LLLLAAWLAE NGRHHECEQL LAWHLFPWSS RFLDVFDHA GHPFYQALGQ LARLTAAQWQ AQLIIPVAVK PLFR
Source	Mammalian cell
Target Names	dmsD
Protein Names	Recommended name: Twin-arginine leader-binding protein DmsD
Expression Region	1-204
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.