



Recombinant Salmonella paratyphi A S-adenosylmethionine synthase (metK)

Product Code	CSB-BP728733SBM
Abbreviation	metK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5PJJ2
Product Type	Recombinant Protein
Immunogen Species	Salmonella paratyphi A (strain ATCC 9150 / SARB42)
Purity	≥85% (SDS-PAGE)
Sequence	MAKHLFTSES VSEGHDPKIA DQISDAVLDA ILQQDPKARV ACETYVKTGM VLVGGEITTS AWVDIEEITR NTVREIGYVH SDMGFDANSC AVLSAIGKQS PDINQGVDRD DPLEQGAGDQ GLMFGYATNE TDVLMPAPIT YAHRLVQRQA EVRKNGTLPW LRPDAKSQVT FQYDDGKIVG IDAVVLSTQH AEDIDQKSLQ EAVMEEIIP ILPSEWLNTS TKFFINPTGR FVIGGPMGDC GLTGRKIIVD TYGGMARHGG GAFSGKDPSK VDRSAAYAAR YVAKNIVAAG LADRCIEQVS YAIGVAEPTS IMVETFGTEK VPAEQLILLV REFFDLRPYG LIQMLDLLHP IYKETAAYGH FGRENFPWEK TDKAQLLRDA AGLK
Source	Baculovirus
Target Names	metK
Protein Names	Recommended name: S-adenosylmethionine synthase Short name= AdoMet synthase EC= 2.5.1.6 Alternative name(s): MAT Methionine adenosyltransferase
Expression Region	1-384
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.