



Recombinant Rat LIX1-like protein (Lix1I)

Product Code	CSB-EP728793RA
Abbreviation	Lix1I
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5PQQ7
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	METMRAQLRQ PGVGVGGRGT LRALRPGVTG APTSAATPPA GPPPAPPPPA PPPPPLLLAG APGLPLPPGA TGSPAVLREA VEAVVRSFAK HTQGYGRAVN VVEALQEFWQ MKQSRGADLK NGALVVYEMV PSNSPPYVCY VTLPGGSCFG SFQFCPTKAE ARRSAAKIAL MNSVFNEHPS RRITDEFIEK SVSEALASFN GNREEADNPN TGIGAFRFML ESNKGKSMLE FQELMTVFQL LHWNGSLKAM RERQCSRQEV LAHYSHRALD DDIRHQMALD WVSREQSVPG ALSRELASTE RELDEARLAG KELRFHKEKK DILMLAAGQL GNMHSSSC
Source	E.coli
Target Names	Lix1I
Protein Names	Recommended name: LIX1-like protein
Expression Region	1-338
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.