



Recombinant Rat PTB domain-containing engulfment adapter protein 1 (Gulp1)

Product Code	CSB-BP728794RA
Abbreviation	Gulp1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5PQS4
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MNRAFSRKKD KTWMHMPEAL SKHYVPYNAK FLGSTEVEQP KGTEVVVRDAV RKLKFARHIK KSEGQKIPKV ELQISYGVK ILEPKTKEVQ HNCQLHRISF CADDKTDKRI FTFICKDSES NKHLCFVFDS EKCAEEITLT IGQAFDLAYR KFLESGGKDV ETRKQIAGMQ KRIQDLETEN MELKNKVQDL ENRLRTTQVS STPAGGVTVT SPSTDIFDMI PFSPISHPSP TSAGNGTQLP PIPSRSEIK RDLFGAEPFD PFNCGAGDFP PDIQSKLDEM QEGFKMGLTL EGTVFCLDPL DSRC
Source	Baculovirus
Target Names	Gulp1
Protein Names	Recommended name: PTB domain-containing engulfment adapter protein 1 Alternative name(s): Cell death protein 6 homolog PTB domain adapter protein CED-6 Protein GULP
Expression Region	1-304
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.