



# Recombinant Rat Inactive serine protease 35 (Prss35)

<b>Product Code</b>	CSB-BP728872RA
<b>Abbreviation</b>	Prss35
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5R212
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	SET EPDFTWHLR IPQVVSEKTI HLASPTFQAD AAVKATVCG IECQEELPAP SLSQLEDFLS YETVFENGTR TLTRVKVQGL VLEPTQNSSI KGARPRRRRQ VYGTDSRFSI LDKRFLTTFP FNTAVKLSTG CSGALVSPNH VLTAHCVHD GKDYVKGSKK LRVGVLKMRN KGGRKRRGS RRSRREAESG GQSPEHPQES TTQRP GKKS RGP RVAQGRP SFQWTRVKST HIPKGWARGE NGDPALDFDY ALLELKRAQK QQYMELGVSP TISKLPGGRI HFSGFDNDRD DQLVYRFCVS SEESNDLLYQ YCDAEAGSTG SGIYLRKKEP GQKNWKRKII AVYSGHQWVD VHGVQKDYNV AVRITPLKYA QICLWIHGNA ANCAYG
<b>Source</b>	Baculovirus
<b>Target Names</b>	Prss35
<b>Protein Names</b>	Recommended name: Inactive serine protease 35
<b>Expression Region</b>	18-406
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.