



# Recombinant Pongo abelii Raftlin-2 (RFTN2)

<b>Product Code</b>	CSB-MP728877PYX
<b>Abbreviation</b>	RFTN2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5R458
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	GCGLRKLED PDDSSPGKIF STLRPQVET KTEFAYEYVL LDFTLQASSN PEVIKINSIL DIVTKVEDYY LKGIVGAIH PVIQVVGQRK HLPASCLYRV VLSRLKLSPK NSAAPSGQRR PRLVIEECPL TSEAQTNDAA KELIEKINFA AKRGMKFVGF ISQPYSYKF CNGTNHDGDI ESMLHVRHSS DENCERSWNEG TLSGQSSSEG IEEELHHEG QYPMEQNGSP SSSKSRKGEA SDNKLYTVFN AFDDDSTSWT YQEGILSMKV TRKGSVISTL DADWLELTTF YYKQGLSLID SFVFWETSKG EHLPKSLEGF FIYEEEGSGV PGSSRKGND AIVVEQWTVIE GCEIKTDYGP LLHTLAIEFGW LLTSVLPTPV LRHDSEGNLA TKQIVFLQRP VMWNSAAQTT DKKASRRIKG EDKNKATSRS IGLDTTTPQP AESRHPPEEC RLSPSRECWT KEGRLAQHNS FSGFSSSDSV LRELDDGQFD QEDGVTQVTC M
<b>Source</b>	Mammalian cell
<b>Target Names</b>	RFTN2
<b>Protein Names</b>	Recommended name: Raftlin-2 Alternative name(s): Raft-linking protein 2
<b>Expression Region</b>	2-501
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.