



Recombinant Pongo abelii Pyruvate dehydrogenase E1 component subunit alpha, somatic form, mitochondrial (PDHA1)

Product Code	CSB-EP728879PYX-B
Abbreviation	PDHA1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5R490
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	>85% (SDS-PAGE)
Sequence	F ANDATFEIKK CDLHRLEEGP PVTTVLTRED GLKYRMMQT VRRMELKADQ LYKQKIIRGF CHLCDGQEAC CVGLEAGINP TDHLITAYRA HGFTFTRGLS VREILAEITG RKGCGAKGKG GSMHMYAKNF YGGNGIVGAQ VPLGAGIALA CKYNGKDEVC LTLYGDGAAN QGQIFGAYNM AALWKLPCIF ICENNRYGMG TSVAAAAST DYYKRGDFIP GLRVDGMDIL CVREATRFAA AYCRSGKGPI LMELQTYRYH GHSMSPGVS YRTREEIQEV RSKSDPIMLL KDRMVNSNLA SVEELKEIDV EVRKEIEDAA QFATADPEPP LEELGYHIYS SDPPFEVRGA NQWIKFKSVS
Source	E.coli
Target Names	PDHA1
Protein Names	Recommended name: Pyruvate dehydrogenase E1 component subunit alpha, somatic form, mitochondrial EC= 1.2.4.1 Alternative name(s): PDHE1-A type I
Expression Region	30-390
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.