



Recombinant Pongo abelii Uncharacterized protein C20orf4 homolog

Product Code	CSB-MP728911PYX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5R5N9
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	>85% (SDS-PAGE)
Sequence	AAMQMDPEL AKRLFFEGAT VVILNMPRGT EFGIDYNSWE VGPKFRGVKT IPPGIHFLHY SSVDKANPKE VGPRMGFFLS LHQRGLTVLR WSTLREEVDL SPAPESEVEA MRANLQELDQ FLGPYPYATL KKWISLTNFI SEATVEKLQP ENRQICAFSD VLPVLSMKHT KDRMGQNLPR CGTECKSYQE GLARLPEMKP RAGTEIRFSE LPTQMFPEGA TPAEITKHSM DLSYALETVL NKQFPSSPQD VLGELQFAFV CFLLGNYEA FEHWKRLNL LCRSEAAAMK HHTLYINLIS ILYHQLGEIP ADFFVDIVSQ DNFLTSTLQV FFSSACSIIV DATLRKKAEK FQAHLTKKFR WDFAAEPEDC APVVVELPEG IEMG
Source	Mammalian cell
Target Names	AAR2
Protein Names	Recommended name: Uncharacterized protein C20orf4 homolog
Expression Region	2-384
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.