



Recombinant Pongo abelii IST1 homolog (IST1)

Product Code	CSB-YP728929PYX
Abbreviation	IST1
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q5R6G8
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	≥85% (SDS-PAGE)
Sequence	<p>MLGSGFKAER LRVNLRRLVIN RLKLLLEK KKT ELAQKARKEI ADYLAAGKDE RARIRVEHII REDYLVEAME ILELYCDLLL ARFGLIQSMK ELDSGLAESV STLIWAAPRL QSEVAELKIV ADQLCAKYSK GYGKLCRTNQ IGTVNDRLMH KLSVEAPPKI LVERYLIEIA KNYNVPYEPD SVVMAEAPPG VETDLIDVGF TDDVKKGGPG RGGGGGFTAP VGGPEGTVPM PMPMPMPAN TPFSYPLPKG PSDFNGLPMG TYQAFPIIHP PQIPATPPSY ESVDDINADK NISSAQIVGP GPKPEASAKL PSRPADNYDN FVLPELPSVP DTLPTASAGA STSASEDIDF DDLSRRFEEL KKKT</p>
Source	Yeast
Target Names	IST1
Protein Names	Recommended name: IST1 homolog
Expression Region	1-364
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>