



Recombinant Pongo abelii Selenoprotein P (SEPP1)

Product Code	CSB-EP728976PYX-B
Abbreviation	SEPP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5R8W9
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	≥85% (SDS-PAGE)
Sequence	E SQDQSSSLCKQ PPAWSIRDQG PMLNSNGSVT VVALLQASUY LCILQASKLE DLRVCLKKEG YSNISHIVVN HQGISSRLKY THLKNKVSEH IPVYQQEENQ TDVWTLNNGS KDDFLIYDRC GRLVYHLGLP FSFLTFPYVE EAVKIAYCEK KCGYCSLTTL KDEDCKSVS LATVDKTVEA PSPHYHHEHH HNHRHQHLGS SELSENQQPG APDAPTHPAP PGLHHHHKHK GQHRQGH PEN RDMPGSEDIQ DLQKKLCRKR CINQLLCKLP KDESELAPRSC CCHCRHLIFE KTGSAITUQC KENLPSLCSU QGLRAEENIT ESCQURLPPP AUQISQQLIP TEASTSURUK NQAKKUEUPS N
Source	E.coli
Target Names	SELENOP
Protein Names	Recommended name: Selenoprotein P Short name= SeP
Expression Region	20-381
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.