



# Recombinant Pongo abelii E3 ubiquitin-protein ligase RING2 (RNF2)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP728992PYX-B   |
| <b>Abbreviation</b>      | RNF2  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q5R9J5  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)   |
| <b>Purity</b>            | ≥85% (SDS-PAGE)   |
| <b>Sequence</b>          | SQAVQTNGT QPLSKTWELS LYELQRTPQE AITDGLEIVV SPRSLHSELM<br>CPICLDMLKN TMTTKECLHR FCADCIITAL RSGNKECPTC RKKLVSKRSL<br>RPDPNFDALI SKIYPSRDEY EAHQERVLAR INKHNNQQAL SHSIEEGLKI<br>QAMNRLQRGK KQQIENGSGA EDNGDSSHCS NASTHSNQEA<br>GPSNKRTKTS DDSGLEPDNN NAAMAIDPVM DGASEIELVF RPHPTLMEKD<br>DSAQTRYIKT SGNATVDHLS KYLAVRLALE ELRSKGESNQ MNLDTASEKQ<br>YTIYIATASG QFTVLNGSFS LELVSEKYWK VNKPMELYA PTKEHK           |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | RNF2  |
| <b>Protein Names</b>     | Recommended name: E3 ubiquitin-protein ligase RING2 EC= 6.3.2.- Alternative name(s): RING finger protein 1B Short name= RING1b RING finger protein 2  |
| <b>Expression Region</b> | 2-336   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Full Length of Mature Protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |