



Recombinant Pongo abelii RNA-binding protein NOB1 (NOB1)

Product Code	CSB-BP729023PYX
Abbreviation	NOB1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5RBB3
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	≥85% (SDS-PAGE)
Sequence	APVEHVVD AGAFLRDAAL QDIGKNIYTI REVVTEIRDK ATRRRRLAVLP YELRFKEPLP EYVRLVTEFS KKTGDYPSLS ATDIQVLALT YQLEAEFVGV SHLKQEPQKV KVSSSIQHPE TPLHISGFHL PSKPKPPQET EKGHPACEPE NLEFSSFMFW RNPLPNIDHE LQELLIDRGE DIPSEEEEE NGFEDRRDDS DDGGGWITP SNIKQIQQEL EQCDVPEDVR VGCVTDFAM QNVLLQMGLH VLAVNGMLIR EARSYILRCH GCFKTTSDMS RVFCSHCGNK TLKKVSVTVS DDGTLHMHFS RNPKVLNPRG LRYSLPTPKG GKYAINPHLT EDQRFQQLRL SRKARQKTNV FAPDYVAGVS PFVENDISSR SATLQVRDSS LGAGRRRLNP NASRRKFVKK R
Source	Baculovirus
Target Names	NOB1
Protein Names	Recommended name: RNA-binding protein NOB1
Expression Region	2-411
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.