



Recombinant Pongo abelii RWD domain-containing protein 3 (RWDD3)

Product Code	CSB-MP729066PYX
Abbreviation	RWDD3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5RDC7
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	≥85% (SDS-PAGE)
Sequence	MAEPVQEELS ALAAIFCGPH EWEVLSRSET DGTVFRIHTK AEGFMDADIP LESVFHLPVN YPSCLPGISI NSEQLTRAQC VTVKEKLLEQ AESLLSEPMV HELVLWIQQN LRHILNQPET GSGYEKCTFS TSTTVDDGLW ITLLHLDHMR AKTKYVKTVE KWASDLRLTG RLMFMGKIIL ILLQGDRNNL KEYLILQKTS KVDVDSSGKK CKEKMISVLF ETKVQTEHKR FLAFEVKEYS ALDELQKEFE TAGLKCLFSE FVLALVK
Source	Mammalian cell
Target Names	RWDD3
Protein Names	Recommended name: RWD domain-containing protein 3
Expression Region	1-267
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.