



Recombinant Danio rerio Sorting nexin-17 (snx17)

Product Code	CSB-MP729124DIL
Abbreviation	snx17
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5RID7
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MHFSIPETEV RSDENGSSYV AYNIHVNGVL HCRVRYSQLL GLHEQIKKEY GNNVVPAPFP KKIFTLTPAE VDQRREQLEK YMQAVRQDPI LGSSEMFNSF LRKAQQETQQ IPTEEVQLEI YLSNGQKVKV NILTSDQTED VLEAVASKLD LPDELVGYS LFLVQERADG SCTYVRKLQE FELPYVSITS LHSSDYRIIL RKS YWDTAYD SDVMDDRVL NLLYAQTVSD IDRGWILVNK EQHRQLKSLQ EKGSKKEFIR LAQTLKYYGY IKFDPCITDF PEKGCHVIVG AGNNELNFHV KLPNEQMKEG SFKVTRMRCW RVTSSQVPVA NGTANPSSSS KCDVKLELAF EYLMKDRQLQ WVTITSQQAI MMSICLQSMV DELMVKKSSG SIKKQMCKKR LNGSLQRSDS QQAVKSPIL DSPDPNREQV VKLSTKLSSV SLRGLSSSNS AGDISGNDFH GNYAFEGIGD DDL
Source	Mammalian cell
Target Names	snx17
Protein Names	Recommended name: Sorting nexin-17
Expression Region	1-473
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.