



Recombinant Rat Protein FAM71F2 (Fam71f2)

Product Code	CSB-EP729127RA-B
Abbreviation	Fam71f2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5RJN7
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MSKIRGLPPA IRDPGPGVEL GVENGLLQQL IHSPEFNLFS DSVVFESTFI QVTKQGNWMD AYERSATIIL GVTSSVPSLP LPNILLMANV TWPQGPSTC STLGAPVITL SRILPLKYVE LQIYDRTQRI LRVRTVTEKI YYLRLHEKHP QAVFRFWIRL VKILQKGLSI TTKDPRIHFT HCLVPKMPNS STETTPESSR PASSQSSETI MLLAAERAGS SVLDLSNRHR FSEDRYTDTK TDNSSNCKTA SPVASLIGMP MRATLIHSLW EQEDSCETLL QAPVASSLGD NFLGP
Source	E.coli
Target Names	Fam71f2
Protein Names	Recommended name: Protein FAM71F2 Alternative name(s): Protein FAM137B
Expression Region	1-295
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.