



Recombinant Mouse TOM1-like protein 2 (Tom1l2)

Product Code	CSB-MP729184MO
Abbreviation	Tom1l2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5SRX1
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MEFLLGNPFS TPVGQCLEKA TDGSLQSEDW TLNMEICDII NETEEGPKDA IRALKKRLSG NRNYREVMLA LTVLETCVKN CGHRFHLLVA NRDFIDSVLV KIISPKNNPP TIVQDKVLAL IQAWADAFRS SPDLTGVVHI YEELKRRGIE FPMADLDALS PIHTPQRSVP EMDPAATIPR SQTQPRRTAG TYSSPPPASY STLQAPALSV TGPITANSEQ IARLRSELDI VRGNTKVMSE MLTEMVPGQE DSSDLELLQE LNRTCRAHQH RIVELISRVS NEEVTEELLH VNDDLNNVFL RYERFERYRS GRSVQNASNG VLSEVTEDEL IDLGGPSPAV VSPMVGSTAP PSSLSSQLAG LDLGTESVSG TLSSLQQCKP QDGFDMFAQT RGNSLAEQRK TVTYEDPQAV GGLASALDNR KQNSEMIPVA QPSVMDDIEV WLRTDLKGDD LEEGVTSEEF DKFLEERAKA AETVPDLPSP PTEAPAPASN TSTRKKPERS DDALFAL
Source	Mammalian cell
Target Names	Tom1l2
Protein Names	Recommended name: TOM1-like protein 2 Alternative name(s): Target of Myb-like protein 2
Expression Region	1-507
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.