



Recombinant Human Uncharacterized protein C9orf50 (C9orf50)

Product Code	CSB-MP729196HU
Abbreviation	C9orf50
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5SZB4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MFWRRLRPGA QDLAPKGLPG DGDFRRSSDP RLPKLTTPPAL RAALGARGSG DWRIPGGGAA WWPEGDAKPG VGVGRLPPRL PALLTATTRA VRKRGLLRSL LPPPLLSAGA SRESAPRQPG PGERERPRRR VAREDPDFLG AFLGELLPSR FREFLHQLQE KCAEEPEPLT SPAPQHQRGV LEHCPGSPRC PNCSFLPDLW GQSSHLQDSL TKISLQQTPI LGPLKGDHSQ FTTVRKANHR PHGAQVPRLK AALTHNPSGE GSRPCRQRCP FRVRFADETL QDITLRYWER RRSVQQSVIV NQKAALPVAS ERVFGSVGKR LESLPKALYP GAKEETLASS SCWDCAGLST KKTQGYLSED TSMNSSLPFC SWKKAAAQRP RSSLRAFLDP HRNLEQESLL PNRVLQSVLK QGCPKGYHLL LASATLQPKD R
Source	Mammalian cell
Target Names	C9orf50
Protein Names	Recommended name: Uncharacterized protein C9orf50
Expression Region	1-431
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.