



Recombinant Human G patch domain-containing protein 4 (GPATCH4)

Product Code	CSB-YP729203HU
Abbreviation	GPATCH4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5T3I0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MNVTPEVKSR GMKFAEEQLL KHGWTQGKGL GRKENGITQA LRVTLKQDTH GVGHDPKEF TNHWWNELFN KTAANLVVET GQDGVQIRSL SKETTRYNHP KPNLLYQKFV KMATLTSGGE KPNKDLESCS DDDNQGSKSP KILTDEMLLQ ACEGRTAHKA ARLGITMKAK LARLEAQEQA FLARLKGQDP GAPQLQSESK PPKKKKKKRR QKEEEEATAS ERNDADEKHP EHAEQNIRKS KKKKRRHQEG KVSDEREGTT KGNEKEDAAG TSGLGELNSR EQTNQSLRKG KKKKRWHHEE EKMGVLEEGG KGKEAAGSVR TEEVESRAYA DPCSRKRKRQ QQEEEDLNLE DRGEETVLGG GTREAESRAC SDGRSRKSKK KRQQHQEEED ILDVRDEKDG GAREAESRAH TGSSSRGKRK RQQHPKKERA GVSTVQKAKK KQKKRD
Source	Yeast
Target Names	GPATCH4
Protein Names	Recommended name: G patch domain-containing protein 4
Expression Region	1-446
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.