



Recombinant Human Putative DNA repair and recombination protein RAD26-like (RAD26L)

Product Code	CSB-BP729209HU
Abbreviation	RAD26L
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5T890
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MQPGSAPPPG RMDPSAPQPR AETSGKDIWH PGERCLAPSP DNGKLCEASI KSITVDENGK SFAVVLYADF QERKIPLKQL QEVKFKDCP RNLIFDDEDL EKPYFPNRKF PSSSVAFKLS DNGDSIPYTI NRYLRDYQRE GTRFLYGHYI HGGGCILGDD MGLGKTVQVI SFLAAVLHKK GTREDIENNM PEFLLRSMKK EPLSSTAKKM FLIVAPLSVL YNWKDELDTW GYFRVTVLHG NRKDNELIRV KQRKCEIALT TYETLRLCLD ELNSLEWSAV IVDEAHRIKN PKARVTEVMK ALKCNVRIGL TGTILQNNMK ELWCVMDWAV PGLLGSGTYF KKQFSDPVEH GQRHTATKRE LATGRKAMQR LAKKMSGWFL RRTKTLIKDQ LPKKEDRMVY CSLTDFQKAV YQTVLETEDV TLILQSSEPC TCRSGQKRRN CCYKTNSHGE TVKTLYLSYL TVLQKVANHV ALLQAASTSK QQETLIKRIC DQVFSRFPDF VQKSKDAAFE TLSDPKYSGK MKVLQQLLNH CRKNRDKVLL FSFSTKLLDV LQQYCMASGL DYRRLDGSTK SEERLKIVKE FNSTQDVNIC LVSTMAGGLG LNFVGANVVV LFDPTWNPAN DLQAIDRAYR IGQCRDVKVL RLISLGTVEE IMYLRQIYKQ QLHCVVVGSE NAKRYFEAVQ GSKEHQGELF GIHNLFKFRS QGSCCLKDIL EV
Source	Baculovirus
Target Names	ERCC6L2
Protein Names	Recommended name: Putative DNA repair and recombination protein RAD26-like EC= 3.6.4.-
Expression Region	1-712
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.