



Recombinant Human Putative uncharacterized protein LOC388900

Product Code	CSB-MP729218HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5TG08
Product Type	Recombinant Protein
Purity	≥85% (SDS-PAGE)
Sequence	MEVINLTTLP MIPVDEHLAV SLVARNTMVK TVRKELENNP PSCLIGSMHQ VNQKIADINL RTEPSANSLA IERFELEKKA LREKTRSSPE DKVKRQRKSQ YSCKGSELRH ARSSVIK RKT ADKNLLAELY QYSNFNSSKP NKLPNGVDFC DMVGNNVRAE RDCLSGKHFC SGRELEKFLS SSSPRAIWL D SFWWIFHERY QPNKELQNNL FDRIAQHYAL LLFRVPKSHS EEALLKRLPS LLSKAVYTSF CCCFPQSWFD THEFKSDICN TMSLWISGGT YPSPQSYDSW DYSELDPERF RREELMLYRR RLTKGREFSL FAGKRAFSQK PAQSRKFYHP QSSSANS PSE KTSSAKQNSE KSLRMQNTAK EHHCQTLVLK KPTQEVKRIS EARECENMFP KKSCAACKSP ELTSNLFNIY GKSPLIVYFL QNYASLQQHG KNVLIVRREK TTSTPDCTPT YTDVISETLC SMKKRKDNLN QLYQHHWTEW NYFDKHLKEL QDNFSREMKN IDPKAADTKK ANHMFIPPSA VNEESPDKKT KGSLQREIEF KGCSNKNHGK GRSVNMGRKE EREREKQKL NISFHSPLSP EELHNLEPGS AYRIRDISAT RTPLRRPSSW EPSLIYQEQV LLKTHTLSME EDKFRC
Source	Mammalian cell
Protein Names	Recommended name: Putative uncharacterized protein LOC388900
Expression Region	1-646
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.