



Recombinant *Oryza sativa* subsp. japonica Thioredoxin-like 2, chloroplastic (Os05g0200100, LOC_Os05g11090)

Product Code	CSB-EP729223OFG-B
Abbreviation	Os05g0200100, LOC_Os05g11090
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5TKD8
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	AH AAVSDKSEQP KWWEKNAGPN MIDIHSTQEF LDALRDAGDR LVIVEFYGTW CGSCRALFPR LCRTAVENPD ILFLKVNFDK NKPMCKRLNV KVLPHYFHFYR GADGQLEAFS CSLAKFQKLK DAIHVHNTAR CSIGPPVGVG DVLDSPEEKP AEASPR
Source	<i>E. coli</i>
Target Names	Os05g0200100
Protein Names	Recommended name: Thioredoxin-like 2, chloroplastic Alternative name(s): Lilium-type thioredoxin 2 OsTrx19
Expression Region	59-216
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.