



# Recombinant Rat Dimethyladenosine transferase 2, mitochondrial (Tfb2m)

<b>Product Code</b>	CSB-EP729242RA-B
<b>Abbreviation</b>	Tfb2m
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5U2T7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	SDLYTQ PLPDCDFEES SSWTHKSRSE PTRHIACKKS ARNLVRDLLE HQNPShQLIL ECNPGPGILT GALLKAGARV VAFESEKMF I PHLESLRKNA DGELQVVHCD FFKIDPRYQE LVRPDVNSHT IFQNLGIKAV PWSAGVPIKV FGILPNKHER RLLWKILFDL YSCESIYRYG RVELNMFISE KEFRKLIATP KRPDLYQVLG VLWQVACEIK FLHMEPWSSF SVHAENGHLE KSKHSESLNL LKQNLVLRM TPRRTLFTEN LSPLNYDMFF HLVKHCFGKR NAPIIHLRS LSTVDPINIL RQIRKRP GDT AAKMYPHDFK RLFETIERSE DSVFKWIYDY CSDDSEL
<b>Source</b>	E.coli
<b>Target Names</b>	Tfb2m
<b>Protein Names</b>	Recommended name: Dimethyladenosine transferase 2, mitochondrial EC= 2.1.1.- Alternative name(s): Mitochondrial 12S rRNA dimethylase 2 Mitochondrial transcription factor B2 Short name= mtTFB2 S-adenosylmethionine-6-N', N'-ade
<b>Expression Region</b>	45-397
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.