



# Recombinant Danio rerio Protein LTV1 homolog (Itv1)

<b>Product Code</b>	CSB-BP729249DIL
<b>Abbreviation</b>	Itv1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5U3J8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Danio rerio (Zebrafish) (Brachydanio rerio)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MPHRKKKSF I NKKNAVSFHL VHRSQKDPLA ADETAPQHVL LPTVRVELEK RKEEQRNFGV FFDDDYDYLQ HLRESTQTAE LISAPRVRRD TRLKPADDEE HMKEEDSSVS EASIKLPSSV FASEFEEVVG LLNKAAPISG PRLDMDPDIV AALDEDFDFD DPENMLEDDF IIKASDVMGG GGGEDDDEWE DTDDDDGEEQ DFDSEAGLSD DEEGGRREFM FADCETKTRF TEYSLTSSVM RRNEQLTLLD DRFDEFYKQF DDDEIGALDN AELEGYIEPD SARLEEVIKD YFIQKEKDCQ KPDQLGPAEL PSVREEDEED EDADEELETE RMVIEPPAER WDCETIISTY SNLYNRPKLI QDPPKPKQIR VSSKSGIPLD VLPKRGLTVK QVERMERIND SDLPRVSTQP RSRDENTEER KARKQAIKSE RKERRTKKKA NKLAFKQEKQ MQEKQMVHLR ANVQGMKLS
<b>Source</b>	Baculovirus
<b>Target Names</b>	Itv1
<b>Protein Names</b>	Recommended name: Protein LTV1 homolog
<b>Expression Region</b>	1-469
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.