



# Recombinant Human RNA-binding protein MEX3C (MEX3C)

<b>Product Code</b>	CSB-EP729263HU
<b>Abbreviation</b>	MEX3C
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5U5Q3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MPSGSSAALA LAAAPAPLPQ PPPPPPPPPP PLPPPSGGPE LEGDGLLLRE RLAALGLDDP SPAEPGAPAL RAPAAAAQGG ARRAAELSPE ERAPPGRPGA PEAAELELEE DEEEGEEAEL DGDLLLLLEEL EEAEEDRSS LLLLSPPAAT ASQTQQIPGG SLGSVLLPAA RFDAREAAAA AAAAGVLYGG DDAQGMMAM LSHAYGPGGC GAAAAALNGE QAALLRRKSV NTTECVVPS SEHVAEIVGR QGCKIKALRA KTNTYIKTPV RGEPIFVVT GRKEDVAMAK REILSAAEHF SMIRASRNKN GPALGGLSCS PNLPGQTTVQ VRVPYRVVGL VVGPKGATIK RIQQQTHTYI VTPSRDKEPV FEVTGMPENV DRAREEIMH IAMRTGNYIE LNEENDFHYN GTDVSFEGGT LGSAWLSSNP VPPSRARMIS NYRNDSSSSL GSGSTDSYFG SNRLADFSPT SPFSTGNFWF GDTLPSVGSE DLAVDSPAFD SLPTSAQTIW TPFEPVNPLS GFGSDPSGNM KTQRRGSQPS TPRLSPTFPE SIEHPLARRV RSDPPSTGNH VGLPIYIPAF SNGTNSYSSS NGGSTSSSPP ESRRKHDCVI CFENEVIAAL VPCGHNLFCM ECANKICEKR TPSCPVCQTA VTQAIQIHS
<b>Source</b>	E.coli
<b>Target Names</b>	MEX3C
<b>Protein Names</b>	Recommended name: RNA-binding protein MEX3C Alternative name(s): RING finger and KH domain-containing protein 2 RING finger protein 194
<b>Expression Region</b>	1-659
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This gene encodes a member of a family of proteins with two K homology (KH) RNA-binding domains and a C-terminal RING-finger domain. The protein interacts with mRNA via the KH domains, and the protein shuttles between the nucleus and cytoplasm. Polymorphisms in this gene may contribute to



hypertension.

---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.