



Recombinant Human RNA-binding protein MEX3C (MEX3C)

Product Code	CSB-MP729263HU
Abbreviation	MEX3C
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5U5Q3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MPSGSSAALA LAAAPAPLPQ PPPPPPPPPP PLPPPSGGPE LEGDGLLLRE RLAALGLDDP SPAEPGAPAL RAPAAAAQGG ARRAAELSPE ERAPPGRPGA PEAAELELEE DEEEGEEAEL DGDLLLLLEEL EEAEEDRSS LLLLSPPAAT ASQTQQIPGG SLGSVLLPAA RFDAREAAAA AAAAGVLYGG DDAQGMMAM LSHAYGPGGC GAAAAALNGE QAALLRRKSV NTTECVVPS SEHVAEIVGR QGCKIKALRA KTNTYIKTPV RGEPIFVVT GRKEDVAMAK REILSAAEHF SMIRASRNKN GPALGGLSCS PNLPGQTTVQ VRVPYRVVGL VVGPKGATIK RIQQQTHTYI VTPSRDKEPV FEVTGMPENV DRAREEIMH IAMRTGNYIE LNEENDFHYN GTDVSFEGGT LGSAWLSSNP VPPSRARMIS NYRNDSSSSL GSGSTDSYFG SNRLADFSPT SPFSTGNFWF GDTLPSVGSE DLAVDSPAFD SLPTSAQTIW TPFEPVNPLS GFGSDPSGNM KTQRRGSQPS TPRLSPTFPE SIEHPLARRV RSDPPSTGNH VGLPIYIPAF SNGTNSYSSS NGGSTSSSPP ESRRKHDCVI CFENEVIAAL VPCGHNLFCM ECANKICEKR TPSCPVCQTA VTQAIQIHS
Source	Mammalian cell
Target Names	MEX3C
Protein Names	Recommended name: RNA-binding protein MEX3C Alternative name(s): RING finger and KH domain-containing protein 2 RING finger protein 194
Expression Region	1-659
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a member of a family of proteins with two K homology (KH) RNA-binding domains and a C-terminal RING-finger domain. The protein interacts with mRNA via the KH domains, and the protein shuttles between the nucleus and cytoplasm. Polymorphisms in this gene may contribute to



hypertension.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.