



Recombinant Acanthamoeba polyphaga mimivirus

Uncharacterized protein R665 (MIMI_R665)

Product Code	CSB-MP729272ADAZ
Abbreviation	MIMI_R665
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5UNS8
Product Type	Recombinant Protein
Immunogen Species	Acanthamoeba polyphaga mimivirus (APMV)
Purity	≥85% (SDS-PAGE)
Sequence	MEHTEPEKTD LKSVDLIEME DTEKLIGTKR KFTDVPVIFR PKKKRVIIIEY EPNTLDIPVS DQESFLKLSA PYITQLYSAP NVTEYTLRDF FVRDKIMPGL IQIQEKMFNI KKTSIRPCYI CKRPMTKVHW FYLYQCNECG NESLKRYRYAS RNLYGKNALV IGGRIKLG YQ IALKLLRAGC RVMITTRSVH NALEFYNQEP DYNEFKSRLF VYSKGFDLGK ARESLPDLVN EITNLF GQQQ LDILIQNAAQ TICFQENIDN KQQEFESLSE AEILD LVKKR GGKRLTFPPV DWRVDFAERY GRVLDRRSIN TWSKNIYETS DEEISVIQS NMIGSILIDK YLIPNMRPSE DTYVIHVHAK EGTFDTHKTM KHTHTNIAKA GLAMLTRCLA GQSENDTEPR SQWPQIHGVN PGWFSIDEYT IYARTRGKIF NPPIDEIDAA SRVVHPIFFQ KKS YCKSWID YKKCHHF
Source	Mammalian cell
Target Names	MIMI_R665
Protein Names	Recommended name: Uncharacterized protein R665
Expression Region	1-467
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.