



# Recombinant Acanthamoeba polyphaga mimivirus Putative ankyrin repeat protein R784 (MIMI\_R784)

<b>Product Code</b>	CSB-EP729315ADAZ
<b>Abbreviation</b>	MIMI_R784
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5UPS1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Acanthamoeba polyphaga mimivirus (APMV)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MINWLHLPDE IWLLIGNFIS DREICCLAFV NTRFLNLLSY NKINQNLVLE CAIEVDNIYI INKLIDKTDV NGLKFVMRWS ALHGSEFQIIK SIIRKDHNNND LLDLHAINWS IESNHLHISK YLIDKYNIVI IPESAFIYGR LDIIDWMFSR QNQNNANKKS TIKKTVDFLK SKFISHDYFY DSMIELACLY GHVHVYQYFI DNQCSVRQKW LYYSCLSGNF DLVLLCKNG CRIFSNRRLI NTAITGGNLD IVRYCLLHISK IDLAQNNFAM KIAIKTGHIG IVRLLVSHGS DIHFDNGECM IIASRGGFAN IVKFFLENKV YMSEKVLKIA AIRGYLDIIK VFIKYASACM SNINSVCKSN GSMVVDNHIN YIVNITRNFN MRKILIWSLI NNRINIVTYL LEIFPKLREI LNQIMYQHPE ISEKIDLDSL YKN
<b>Source</b>	E.coli
<b>Target Names</b>	MIMI_R784
<b>Protein Names</b>	Recommended name: Putative ankyrin repeat protein R784
<b>Expression Region</b>	1-433
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.