



Recombinant Acanthamoeba polyphaga mimivirus Uncharacterized protein R255 (MIMI_R255)

Product Code	CSB-MP729317ADAZ
Abbreviation	MIMI_R255
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5UPT9
Product Type	Recombinant Protein
Immunogen Species	Acanthamoeba polyphaga mimivirus (APMV)
Purity	≥85% (SDS-PAGE)
Sequence	MSDTIYNPKI ITIEQFESFS ERKVKPIDGL NIPNIFYPPH LVNKKRPTEA DNSLTINIRQ LFNSLSDTNI NKVKDQLKET IITKAKNENM IEEIAKEILS NFVISEKNIR NYMVLLNAVS GTCVLLADSE TKKTSPTIGK YFIDNCKTDI FHHISESTVR ELANMDLDDS DQLDLYNRKR ETVINLILTL CFLYGQRNTD LIRLTASQLY PLINTIMNTY DSLQAKMKIL GNPYDGGDDCE DEEEYEILSK MCTIYAEQLY TFINKEVSSF LLDDTVVKGM LMKNLVERFK NTIVPTLTEA YLRKCSICIQ YN
Source	Mammalian cell
Target Names	MIMI_R255
Protein Names	Recommended name: Uncharacterized protein R255
Expression Region	1-312
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.