



Recombinant Acanthamoeba polyphaga mimivirus Putative ankyrin repeat protein R787 (MIMI_R787)

Product Code	CSB-MP729327ADAZ
Abbreviation	MIMI_R787
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5UQ08
Product Type	Recombinant Protein
Immunogen Species	Acanthamoeba polyphaga mimivirus (APMV)
Purity	≥85% (SDS-PAGE)
Sequence	MDSMLIVNPF DALDQEFYHS LPPELWVVLL NTSIKSSINL AFTCKQFFQL CFLKSSSTNV IKLAVRNGCI HILKYIDKIK LIDTSHSKNT RFELSDIQNY LRKACKYGHL SVVKYLVGRG ADIKAEKSYL TGPITNGHLE VVKYMVSMGV DFRYDDDAML REAALKGQFE MVKYLCEIGS DVSSKNNLTL VWAIIGGHID IVDYLSKGA DYRQINNPLA WACHHGHFNI VEYFVDKNVD IKADDNLALR QAAENGHLSI VKYLVDLGAD INAKNSCAIR WSSRQGHNLV VQFLINQKAD VNAKDNYAIN HACKKGYTDI IKCLVESGAE FKHNNNYCLR IAAGKGHMSI VRYLVEKGS D YTDLDNYAIK SALCNCKFIV AKYLINCGVD LRDNNAYILR YCLKKSSIPT KLRIVKFALD TYHDVNLITN CFASITK
Source	Mammalian cell
Target Names	MIMI_R787
Protein Names	Recommended name: Putative ankyrin repeat protein R787
Expression Region	1-437
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.