



Recombinant Acanthamoeba polyphaga mimivirus Putative FNIP repeat-containing protein R636 (MIMI_R636)

Product Code	CSB-YP729378ADAZ
Abbreviation	MIMI_R636
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5UR85
Product Type	Recombinant Protein
Immunogen Species	Acanthamoeba polyphaga mimivirus (APMV)
Purity	≥85% (SDS-PAGE)
Sequence	MIDHFIMNKL PLELVYYLTN WLKDYDKRML LSTCKNFLSL QSCVYYEDIY YYEEIHKLSY HQMFRKVELE LCKFKTPIPK IVTNLEFDVY FNEPLINRRI KIRNKWVKIR KIIPSTVKYL DMGYSFNKSI RGALSYGLTT LIFGPEFNQP IDNYIPQTVT YLEFSVYFDQ PVVGYLPTL THLIFGTDFN QPIKGALPDT LEYLYFGWAF NQPIDFALPP NLKCLKFGHC FDQPIKGFLP LGLEKLIFEA EYHQPMDLAI PESVEYLKLG NPGSNSLENC LPINLKTUVF GENFNENLNT LSLDKFRELE KIEFDGFFDQ KISSRITIDN VIKSNNFIET FQNFYGMIKY YFNDFINPNR QYLPENVKHI IISTYLYDKI CSNISNDVFI EVYN
Source	Yeast
Target Names	MIMI_R636
Protein Names	Recommended name: Putative FNIP repeat-containing protein R636
Expression Region	1-394
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.