



Recombinant Subtilisin-like protease 1 (SUB1)

Product Code	CSB-MP729415TBAV
Abbreviation	SUB1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5VJ77
Product Type	Recombinant Protein
Immunogen Species	Trichophyton verrucosum (Cattle ringworm fungus)
Purity	≥85% (SDS-PAGE)
Sequence	GKVE QSNVPSWGLA RISNPQPGAG SYIYDSSAGE GITVYSVDTG VDVNHEDFEG RAIWGSNQVN DGDDRDRGSGH GTHTSGTMVG KEFGIAKKAK LVAVKVLGND GSGPTSGIVA GINWSVEHAR QNGGTTKAVM NMSLGGSSSS ALNRAAAQAV EQGMFLSVAA GNDNQDAQSS SPASEPSVCT VGSSAEDDSR SSFSNWGPAI DLFAPGSNII SARPGGGSQS MSGTSMAAPH VAGLAAYLMA LEGISGGAVC DRLKELGTSS ITDAGPGTPT NVLINNGGA
Source	Mammalian cell
Target Names	SUB1
Protein Names	Recommended name: Subtilisin-like protease 1 EC= 3.4.21.-
Expression Region	117-399
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.