



Recombinant Human KH domain-containing, RNA-binding, signal transduction-associated protein 2 (KHDRBS2)

Product Code	CSB-BP730139HU
Abbreviation	KHDRBS2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5VWX1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEEEKYLPEL MAEKDSLDPV FVHASRLLA EIEKFQGS DG KKEDEEKYL DVISNKNIKL SERVLIPVKQ YPKFNFV GKL LGPRGNSLKR LQEETGAKMS ILGKGS MRDK AKEEELRKSG EAKYAHLSDE LHV LIEVFAP PGEAYSRMSH ALEEIKKFLV PDYNDEIRQE QLRELSY LNG SEDSGRGRGI RGRGIRIAPT APSRGRGGAI PPPPPPGRGV LTPRGSTVTR GALPVPPVAR GVPTPRARGA PTVPGYRAPP PPAHEAYEEY GYDDGYGGEY DDQTYETYDN SYATQTQSVP EYYDYGHGVS EDAYDSYAPE EWATTRSSLK APPQRSARGG YREHPYGRY
Source	Baculovirus
Target Names	KHDRBS2
Protein Names	Recommended name: KH domain-containing, RNA-binding, signal transduction-associated protein 2 Alternative name(s): Sam68-like mammalian protein 1 Short name= SLM-1 Short name= hSLM-1
Expression Region	1-349
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.