



Recombinant Human Doublesex- and mab-3-related transcription factor A1 (DMRTA1)

Product Code	CSB-BP730143HU
Abbreviation	DMRTA1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5VZB9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MERSQCGSRD RGVSGRPHLA PGLVVAAPPP PSPALPVPSG MQVPPAFLRP PSLFLRAAAA AAAAAAATSG SGGCPPAPGL ESGVGAVGCG YPRTPKCARC RNHGVVSALK GHKRFCRWRD CACA KCTLIA ERQRMMAAQV ALRRQQAQEE SEARGLQRLL CSGLSWPPGG RASGGGGRAE NPQSTGGPAA GAALGLGALR QASGSATPAF EVFQQDYPEE KQE QKESKCE SCQNGQEELI SKSHQLYLGS SSR SNGVIGK QSIGSSISEY SNKPDSILSP HPGEQSGGEE SPRSLSSSDL ESGNESEWVK DLTATKASLP TVSSRPRDPL DILT KIFPNY RRSRLEGILR FCKGDVVQAI EQVLNGKEHK PDNRNLANSE ELENTAFQRA SSFSLAGIGF GTLGNKSAFS PLQTT SASYG GDSSLYGVNP RVGISPLRLA YSSAGRGLSG FMSPYLTPGL VPTLPFRPAL DYAFSGMIRD SSYLSSKDSI TCGRLYFRPN QDNP
Source	Baculovirus
Target Names	DMRTA1
Protein Names	Recommended name: Doublesex- and mab-3-related transcription factor A1
Expression Region	1-504
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.