



Recombinant Bacillus clausii UPF0403 protein ABC2045 (ABC2045)

Product Code	CSB-MP730184BAAC
Abbreviation	ABC2045
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5WGC5
Product Type	Recombinant Protein
Immunogen Species	Bacillus clausii (strain KSM-K16)
Purity	≥85% (SDS-PAGE)
Sequence	MSMAYEEYMR QVVLPMRQEL TDAGFKELTT AEAVEDYIAD AEGTTLVFIN SVCGCAAGLA RPAAIYSLQH DHTPDHLVTV FAGQDKEATA AMRSKFPDYP PSSPSMALLK GTEVVHFIPR EDIEGAEPEA IVRNLALAYN EHCQK
Source	Mammalian cell
Target Names	ABC2045
Protein Names	Recommended name: UPF0403 protein ABC2045
Expression Region	1-145
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.