



Recombinant Dog Lymphotoxin-alpha (LTA)

Product Code	CSB-MP730213DO
Abbreviation	LTA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5WR07
Product Type	Recombinant Protein
Immunogen Species	Canis lupus familiaris (Dog) (Canis familiaris)
Purity	>85% (SDS-PAGE)
Sequence	FPGVGIS PSAARTAYQH PQKHFIQGT LKPA AHLIGDP SIQNSLRWRA NTDRAFLRHG FLSLNNSLLV PTSGLYFVYS QVVFSGEGCF PKATPTPLYL AHEVQLFSSQ YPFHVPLLSA QKSVCPGPQG PWVRSVYQGA VLLTQGDQL STHTDGISHL LLSPSSVFFG AFAL
Source	Mammalian cell
Target Names	LTA
Protein Names	Recommended name: Lymphotoxin-alpha Short name= LT-alpha Alternative name(s): TNF-beta Tumor necrosis factor ligand superfamily member 1
Expression Region	34-204
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.