



Recombinant Legionella pneumophila UPF0133 protein lpp2803 (lpp2803)

Product Code	CSB-MP730267LAAB
Abbreviation	lpp2803
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5X1E0
Product Type	Recombinant Protein
Immunogen Species	Legionella pneumophila (strain Paris)
Purity	≥85% (SDS-PAGE)
Sequence	MDINQNLGNL MKEAQKMQQR MQEAQQQLSQ LVVSGESGGG MVTIKMNGRH DVTEVKIKPT LMDEDIEMLE DLIAAAVNDA VRKIEKASKE KISQLTAGLN IPTDLMGGKE GE
Source	Mammalian cell
Target Names	lpp2803
Protein Names	Recommended name: UPF0133 protein lpp2803
Expression Region	1-112
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.