



Recombinant Metal ABC transporter substrate-binding lipoprotein (mtsA)

Product Code	CSB-YP730342SMZ
Abbreviation	mtsA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5XDI6
Product Type	Recombinant Protein
Immunogen Species	Streptococcus pyogenes serotype M6 (strain ATCC BAA-946 / MGAS10394)
Purity	≥85% (SDS-PAGE)
Sequence	CSSTGKTAK SDKLKVVATN SIIADMTKAI AGDKIDLHSI VPIGQDPHEY EPLPEDVEKT SNADVIFYNG INLEDGGQAW FTKLVKNAQK TKNKDYFAVS DGIDVIILEG ASEK GKEDPH AWLNLENGII YSKNIAKQLI AKDPKNKETY EKNLKAYVAK LEKLDKEAKS KFDAIAENKK LVTSEGCFK YFSKAYGVPS AYIWEINTEE EGTPDQISSL IEKLVKIPKPS ALFVESSVDR RPMETVSKDS GIPIYSEIFT DSIAKKGKPG DSYAMMKWN LDKISEGLAK
Source	Yeast
Target Names	mtsA
Protein Names	Recommended name: Metal ABC transporter substrate-binding lipoprotein
Expression Region	21-310
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.