



Recombinant V-type ATP synthase beta chain (atpB)

Product Code	CSB-YP730350SMZ
Abbreviation	atpB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5XE49
Product Type	Recombinant Protein
Immunogen Species	Streptococcus pyogenes serotype M6 (strain ATCC BAA-946 / MGAS10394)
Purity	≥85% (SDS-PAGE)
Sequence	MSVLKEYRTV SEVVGPLMIV DQVAGVHYNE LVDITLHNGE RRRKGQVLEVQ GDKAMVQLFE GSTGINLAKT KVRFTGHPL E LAVSEDMVGR IFDGMGQPID GGPELIPEKY LDIDGQAINP VARDYPDEFI QTGISAIDHL NTLVRGQKLP VFSGSGLPHN ELAAQIARQA TVLNSDDNFA VVFAAMGITF EEAFFMNDL RETGAIDRSV LFINLANDPA IERIATPRIA LTTAEYLA YE KGMHVLVIMT DMTNYCEALR EVSAARREVP GRRGYPGYLY TNLSTLYERA GRLIGKKGSV TQIPILTMPE DDITHPIPDL TGYITEGQII LSQELYKNGF RPPINVLP SL SRLKDKGSGE GKTRQDHAAT MNQLFAAYA Q GKQAKELAVV LGESALSETD KLYVAFTNRF EEEYINQGFY TNRSIEESLD LGWELLSILP RTELKRIKDD MLDRYLPKAD TTMTKV FVAN D
Source	Yeast
Target Names	atpB
Protein Names	Recommended name: V-type ATP synthase beta chain Alternative name(s): V-ATPase subunit B
Expression Region	1-471
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.