



Recombinant *Xenopus laevis* RNA-binding protein with serine-rich domain 1-A (rnps1-a)

Product Code	CSB-BP730359XBE
Abbreviation	rnps1-a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5XG24
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	≥85% (SDS-PAGE)
Sequence	MAPSPSKRKE RSEDRAKERG KEKAPGKEVT EKDRGRDKAK KRRSGSSGSS SSSHSRSSSS SSSSSGSSSG SSSGSSSSAS SRSGSSSSSR SSSSSSSSGS PPSRRRHDN RRRSRSKSKQ PKRDEKERKR RSPSPRPTKV HIGRLTRNVT KDHILEIFST YGKIKMIDMP VDRYHPHLSK GYAYVEFEAP EEAEKALKHM DGGQIDGQEI TASAVLTPWP MRPMRRFSP PRRMLPPPPM WRRSPPRMRR RSRSPRRRSP VRRRSRSPAR RRHRSRSSSN SSR
Source	Baculovirus
Target Names	rnps1-a
Protein Names	Recommended name: RNA-binding protein with serine-rich domain 1-A
Expression Region	1-283
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.