



# Recombinant Mouse Protein FAM169A (Fam169a)

<b>Product Code</b>	CSB-EP730360MO-B
<b>Abbreviation</b>	Fam169a
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5XG69
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MAFPVDLLDN CTHEELNSS EDYLSSLRCG DPEHPECFSS LNITIPVSLN NVGFVPLYGG NQTQKILALF APEDSLTAVA LYLVGQWWAI DDIVKTSEPS REGLKQVSTL GERVVLYVLN RIIYRKQEME RNEIPFLCHS STDYAKILWK KGEAVGFYSV KPTGSLCASF LTQNYQLPVL DTMFIRKKYR GKDLGLHMLE DFVDSFTEDA LGLRYPLSSL MYTASKQYFE KYPGDHELLW EVEGVGHWHQ RVPVTRALQR EAIKATDVSQ YEATRPVSGE YGLAAVPEHE PGLDDTQSSE LQIHSLKDAF ASTSEGPEKT PVSTRTRSSH LKRPKIGKHF QDSEFSSSQG EDENVAKTSP TASVNKIEYA ARTSESSEEF LEEPEQGI DFEDESGDKD AQALETQPR LQKQDGDKDS ALEPVNGEVM DAALKPSLTT EDEDSTSEGL EEDLKVPPFN SSGEPGNPVP LVAESSKVPE ATLAKTSPDT DSEMLIDQSP SDDKGHTEEN LSPVSKKTL LGSSDNVATV SNIKSDGNF PNSVPEFPE EPVSNLSPN TTSSVEDQGE EGAPEAQEPS ATQSSLIEVE LEDAPFPQNA GQKSQSEEQS EASSEHLEQF TQSAEKAVDS SSEEIEVEVP VVDRRNLRRK AKGHKGP GKK KAKLT
<b>Source</b>	E.coli
<b>Target Names</b>	Fam169a
<b>Protein Names</b>	Recommended name: Protein FAM169A
<b>Expression Region</b>	1-665
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.