



Recombinant Mouse Protein FAM169A (Fam169a)

Product Code	CSB-EP730360MO
Abbreviation	Fam169a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5XG69
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MAFPVDLLDN CTHEELNSS EDYLSSLRCG DPEHPECFSS LNITIPVSLS NVGFVPLYGG NQTQKILALF APEDSLTAVA LYLVGQWWAI DDIVKTSEPS REGLKQVSTL GERVVLYVLN RIIYRKQEME RNEIPFLCHS STDYAKILWK KGEAVGFYSV KPTGSLCASF LTQNYQLPVL DTMFIRKKYR GKDLGLHMLE DFVDSFTEDA LGLRYPLSSL MYTASKQYFE KYPGDHELLW EVEGVGHWHQ RVPVTRALQR EAIKATDVSQ YEATRPSVSGE YGLAAVPEHE PGLDDTQSSE LQIHSLKDAF ASTSEGPEKT PVSTRTRSSH LKRPKIGKHF QDSEFSSSQG EDENVAKTSP TASVNKIEYA ARTSESSEEF LEEPEQGI DFEDESGDKD AQPALETQPR LQKQDGDKDS ALEPVNGEVM DAALKPSLTT EDEDSTSEGL EEDLKVPPFN SSGEPGNPVP LVAESSKVPE ATLAKTSPDT DSEMLIDQSP SDDKGHTEN LSPVSKKKTL LGSSDNVATV SNIEKSDGNF PNSVPEFPE EPVSQNLSPN TTSSVEDQGE EGAPEAQEPS ATQSSLIEVE LEDAPFPQNA GQKSQSEEQS EASSEHLEQF TQSAEKAVDS SSEEIEVEVP VVDRRLRRK AKGHKGP GKK KAKLT
Source	E.coli
Target Names	Fam169a
Protein Names	Recommended name: Protein FAM169A
Expression Region	1-665
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.