



# Recombinant Rat Protein FAM55D (Fam55d)

<b>Product Code</b>	CSB-BP730380RA
<b>Abbreviation</b>	Fam55d
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5XI89
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	RNPV KLWAVFKLPA SFNQWDLIMK SSCPKVPLNP SVSPTETELR IREILEKLNK QIPPRPFAHL NNTTSATHSI ATILNPQDTY CVGDQLDILV EARDHLRNRK GYGGDFLRAR MSSPALKAGA SGKVTDFNNG TYLVSFTLFW EGQVLSILL MHPSEGVSAL WRARNQGYDR IIFSGHFVSG ASQVHTDCAL VLNSSVELCQ YLDAQDQEAF YCVKPPNVPC AAITHMHSKN KDISYLSQQE RSLFERSNIA VEIMGKSNVI SVSKCNKAVP VKKKCKFGMA SAIPTGHVWK NTWNPASCSL APIKMKDCLR GKLVHLMGDS TMRQWMEYFK SKINTLRPVD LHETGRLQHQ LAVDLDEKIN IQWQKHGFPL IGSLVYSVKE IENTARIIDR IGGEKNTVIV FSLGQHFRPF PIDVFIRRAL SVHRALQRLR LRSPDTLVVL KTENTRELNN DMERFSDFHG YTQYLALKNI FQDLRVGVID AWDMTVAYGT NDVHPPEEVV RSEINIFLNY IC
<b>Source</b>	Baculovirus
<b>Target Names</b>	Nxpe4
<b>Protein Names</b>	Recommended name: Protein FAM55D
<b>Expression Region</b>	27-542
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.