



Recombinant Rat Coiled-coil domain-containing protein 86 (Ccdc86)

Product Code	CSB-YP730383RA
Abbreviation	Ccdc86
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5XIB5
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	MGTPLRRSRR LEGLNPLSLE NLPDPEVSRA KRALVDFKSN PEETRELESP RVPLDLVSP QPQPETSPES PCPKQDAGFG SPQRQPEPHP GSLQLHQDLG LDSPAGQTES NPESPQREQS SKLSPTQDSE VAHAKEEVIP RSPEPCPGQQ APGPEPSQPA QELAFQAPSS PERQLEPSKL PPAGESVTGS LDLKRVIAS PQAPASKKLLK EELPVIKPKGK PKSGRVWKDR SKKRF SQMVQ DKPLRTSWQR KMKERQERKL AKDFARHLEE EKQRRRQEKK ERRAENLRRR LENERKAEIV QVIRNPAKLLK KAKKKQLRSI QKRDTLALLQ KQPPQRPVAK V
Source	Yeast
Target Names	Ccdc86
Protein Names	Recommended name: Coiled-coil domain-containing protein 86 Alternative name(s): Cytokine-induced protein with coiled-coil domain
Expression Region	1-341
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.