



Recombinant Human F-box only protein 31 (FBXO31)

Product Code	CSB-BP730405HU
Abbreviation	FBXO31
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5XUX0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAVCARLCGV GPSRGCRRRQ QRRGPAETAA ADSEPDTDPE EERIEASAGV GGGLCAGPSP PPRCSLLEL PELLVEIFA SLPGTDLPSL AQVCTKFRRI LHTDTIWRRR CREEYGVCEN LRKLEITGVS CRDVYAKLLH RYRHILGLWQ PDIGPYGGLL NVVVDGLFII GWMYLPPHDP HVDDPMRFPK LFRIHLMERK AATVECMYGH KGP HHGHIQI VKKDEFSTKC NQTDHHRMSG GRQEEFRTWL REEWGRTLED IFHEHMQELI LMKFIYTSQY DNCLTYRRIY LPPSRPDDLI KPGLFKGTYG SHGLEIVMLS FHGRRARGTK ITGDPNIPAG QQTVEIDLRH RIQLPDLENQ RNFNELSRIV LEVRERVRQE QQEGGHEAGE GRGRQGPRES QPSPAQPRAE APSKGPDGTP GEDGGEPGDA VAAAEQPAQC GGGQPFVLPV GVSSRNEDYP RTCRMCFYGT GLIAGHGFTS PERTPGVFIL FDEDRFVFW LELKSFSLYS RVQATFRNAD APSPQAFDEM LKNIQSLTS
Source	Baculovirus
Target Names	FBXO31
Protein Names	Recommended name: F-box only protein 31
Expression Region	1-539
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.