



Recombinant Nocardia farcinica UPF0345 protein NFA_15710 (NFA_15710)

Product Code	CSB-MP730453NAAA
Abbreviation	NFA_15710
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5YZH4
Product Type	Recombinant Protein
Immunogen Species	Nocardia farcinica (strain IFM 10152)
Purity	≥85% (SDS-PAGE)
Sequence	MSKFENVSVK KKANVYFDGA CVSHSIEFPD GTAKSVGVIL PATLTFGTTA PEVMELVEGH CRVTLPGADA PVTFRRGGESF EVPADSEFTI EVLETVHYVC HYG
Source	Mammalian cell
Target Names	ppnP
Protein Names	Recommended name: UPF0345 protein NFA_15710
Expression Region	1-103
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.