



# Recombinant *Oryza sativa* subsp. japonica IAA-amino acid hydrolase ILR1-like 6 (ILL6)

<b>Product Code</b>	CSB-MP730468OFG
<b>Abbreviation</b>	ILL6
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5Z678
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. japonica (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	TFLEED VILGTVEEAK VARLGGGGGG GSKGANASTR RADNTCAGVG VGVGGGGGGG GGGGGGGRGR FYLGWKEEIA GMAGRPETAA WLRVRRRIH ERPELAYEEV ETSRLVRDEL DAMGVGFRHP VARTGVVANI GTGRPPVVAL RADMDALPIQ EAVEWEHKSK NPGKMHACGH DAHVAMLLGA AKILKAREHH LRGTVRLLFQ PAESGAGAK RMIEGGALED VEAIFAVHVS HQHPTSVIGS RTGPLLAGCG FFKAVIHGGR RSGDAVLA STIISLQSIV SREADPLDSQ VVSVAMVNGS DHPAATARAA AEEEEEEFVL GGTFRAFSNA SFYQVRRRIE EVITAQARVH GCEAAVDFFE NQSFYPPTVN DARMYAHVKA VAGELLGAGS YRDVPPMMGA EDFSFYSQVV PAGFYYIGVR NETLGSVHTG HSPYFMIDED VLPTGAAFHA AIAERYLANH SPSSSSSSSDS DDPDVELEAS
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ILL6
<b>Protein Names</b>	Recommended name: IAA-amino acid hydrolase ILR1-like 6 EC= 3.5.1.-
<b>Expression Region</b>	25-510
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.