



Recombinant Chicken Methyltransferase-like protein 2 (METTL2)

Product Code	CSB-MP730476CH
Abbreviation	METTL2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZHP8
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	MAAPSEAPER RRPFGRRFLT DPTRLFQHNA WDNVEWSEEQ EATAKSKVQE NSSQLLPQDK QEEYEVNAKR YWDDFYKIHE NGFFKDRHWL FTEFPELAPN RNPSQNEDSL CEFSCKEVSK NEGLGSCENG HCTLENRAEN QLNLLKSSPR FCTEELAPQK LKQSYEDYPG SSASYRILEV GCGAGNTVFP ILQTNNDPGL FVYCCDFSTT AVDLVQSNVE YDSSRCFAFV HDLCNDQSPF PMPDESLDIV ILIFVLSAIL PEKMQCVINK LSRLLKPGGM ILLRDYGRYD LAQLRFKKGQ CLSANFYVRG DGTRVYFFTQ DELDDLFTRA GLQKIQNLVD RRLQVNRGKQ MTMYRVWIQC KYQKPAGPQL
Source	Mammalian cell
Target Names	METTL2
Protein Names	Recommended name: Methyltransferase-like protein 2 EC= 2.1.1.-
Expression Region	1-370
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.